



# -150°C Ultra-Low Temperature V.I.P.™ Freezer

## MDF-C2156VAN

- New V.I.P. PLUS™ technology maximises storage capacity.
- Comfortable operation via LCD control panel
- 1. Chamber temperature recorded for up to 1 month
- 2. Warning messages displayed in LCD panel
- 3. Records history of door opening/closing (more than 2 minutes)
- Can hold up to 14,000 tubes (2ml)
- Product conforms to RoHS (European Restriction of Hazardous Substances directives)







Conventional Insulation

SANYO MDF-C2156VAN with V.I.P. Insulation

100%
184L
6.5cu.ft.

100pcs 2" boxes

140pcs 2" boxes

V.I.P.PLUS Evolution series -150°C



#### **Features**

#### Evolutionary Design

- SANYO V.I.P. PLUS<sup>™</sup> Evolution Series -150°C ultra-low temperature freezers use patented revolutionary vacuum insulation cabinet construction that reduces the wall thickness from 175mm to 135mm (6.9" to 5.3") and achieves up to 25% more storage capacity than a conventionally insulated freezer without increasing the footprint.
- The newly developed refrigeration system and the freezer structure offer a quiet experimental environment.
- Achieves 50% higher energy conservation than previous model. (For model rated 220V, 50Hz)

#### Reliability

- Status Alert monitors ambient and system conditions continuously and notifies of any abnormalities before a problem happens.
- Equipped with LN<sub>2</sub> backup system as standard to prepare for any contingency.

#### Ideal replacement for LN2 storage


Compared to liquid nitrogen systems, SANYO MDF-C2156VAN offers the following advantages:

- Lower running cost
- Easy maintenance
- No occurrence of contamination


#### **Performance Data**

#### Pull-down characteristics Pull-up characteristics during power failure


Central air temp.in cabinet Ambient temp. 30°C, No-load



#### LN<sub>2</sub> back up system operation characteristics



#### **Innovative Control Panel**



#### Setting display

T+00 \$+11145		_	# F T W
Temperature	-	1-150.0	0.4
migh bierm		+ 4" 8 - 5 8	
ter Atres	48.4	4874-88	
Bing Bert	10000	10.000	tinted a
Seer Serey	E 401 4	1 4 - 10	00.00
Attent Total		48.000	1.041


#### Access control

To prevent unauthorized access to freezer settings, the key lock can only be cancelled by inputting the password.

#### Data logging function

Recording of internal temperature for 1 month or 1 week period

Recording of history of door opening/closing (more than 2 minutes)




### **Specifications**

Model	MDF-C2156VAN	
Temperature range	-125°C to -150°C (1°C increments)	
Maximum cooling performance	-150°C (Ambient Temp. 30°C)	
Exterior dimensions (W x D x H)	1730 x 765 x 1010 mm (68.1" x 30.1" x 39.8")	
Interior dimensions (W x D x H)	760 x 495 x 615 mm (29.9" x 19.5" x 24.2")	
Net weight	Approx. 319 kg	
Effective capacity	231L (8.2 cu.ft.)	
Access port	40 mm diameter, 1 location	
Compressor	Hermetic type, Output: 1100W (high stage side),	
	1100W (low stage side)	
Refrigerant	HFC	
Alarm	High/Low temperature, Power failure, Filter check, Self diagnosis,	
	Door check, Part replacement notification	
Remote alarm contact	Allowable contact capacity, DC 30V, 2A	
Accessories	1 set of keys, 1 scraper, Connecting tube for LN <sub>2</sub> backup system	

<sup>\*</sup>Appearance and specifications are subject to change without notice

[Unit:mm (inch)]

#### **Dimensions**



#### **Options**



Strip chart recorder MTR-155H

Recorder mounting bracket (for MTR-155H) MDF-S30150 Inventory rack
IR-210C
Max capacity: 14pcs

Chart paper **RP-155** 

Ink pen DF-38FP

SANYO DAQ (Data Acquisition) system

MTR-480/ MTR-2000

Further details regarding this product can be accessed at http://www.sanyo.co.jp/cmg/biomedical

We design and manufacture products in accordance with ISO9001 & 14001 requirements. RoHS (European Restricton of Hazardous Substances directives) compliant

Distributed by:



<sup>\*</sup>Cooling performance is indicated by the temperature reached at the center of the freezer (at ambient temperature of 30°C with no load). In order to use the freezer at a stable temperature for a long time, it is recommended that the temperature be set to at least 5°C higher than the indicated lowest temperature. In addition, depending on the usage conditions, it may not be possible to reach the indicated lowest temperature.