

Nature Sets the Standards

Plant Growth Cabinets and Plant Growth Walk-In Chambers



Bio Line, the system which has been developed "by biologists for biologists". It was created in cooperation with scientists from many different fields of research. Biological aspects, research diversity along with easy operation and cleansing as well as maintenance friendliness have been considered during the development of this line. Thus, *Bio Line* chambers are especially adapted to modern experimental practice in the field of natural science.

Modularity - the concept of *Bio Line*. The system consists of three types of modules: a research room, an electro and machine module and a variety of functional modules. The research room is available in three sizes. The functional modules, e.g. different irradiation intensities, temperature, air flow and height modules adapt the research room to suit not only your experiments but also your plants. Should your research requirements change then *Bio Line* adapts itself to the new requirements by exchanging or retrofitting individual functional modules.

Bio Line offers natural horizontal air flow. An innovative fine pored airflow fiber made of polypropylene ensures not only extremely low air velocity rates in the research room but also high climate constancy.

The low air velocity rate in the *Bio Line* research rooms supports a more natural evaporation rate of the stomata without water stress.

Thanks to its special fibrous structure, the fiber is microbiologically and physiologically safe and very easy to clean.




Fig. 1: Machine module - good access from the front



Fig. 2: *Bio Line* cabinets VB 1014 and VB 1514




Fig. 3: Humidity range
1 = range inapplicable with irradiation

Bio Line - the space saving concept. Access, operation and maintenance from the front. To facilitate service work, the electro and machine module and the connections are accessible via large doors. Power and water connections are easy to connect. Several systems may be installed in one row, shoulder to shoulder, without a gap between them. An innovative humidification system creates germ-free vapour of approx. +150 °C, thanks to shock heated stream. Protein molecules on contaminants denature immediately. This low-maintenance system is operated with demineralized water. A highly efficient separate dehumidification exchanger provides a constant climate.

Excellent accessibility to the research room is provided thanks to large doors with observation windows. Door and window flaps close light-tight. The *Bio Line* research room may be easily cleaned from chemical and micro-biological contamination.

Bio Line does not use materials which contain substances proven to be phytotoxic. The polymeric materials assigned in the research room have been gas-chromatographically analyzed for volatile organic emissions.

The research room is made of insulated elements and is sealed vapour-tight. The insulation material ensures excellent heat insulation which contributes to high constancy of climate conditions and to the reduction of operating costs.

Highly reflecting surfaces create a natural scattering of rays and thus extremely uniform irradiation is achieved over the entire research room.

Bio Line uses the chlorine-free refrigerant. Ozone depletion potential is thus reduced to zero. The insulating material is 50 mm pure mineral fiber which is free from CFC and asbestos.




Fig. 4: Scheme - modularity




Fig. 5: Irradiation module 1




Fig. 6: Airflow-fiber

Technical Data *Bio Line* Cabinets

Type	Fig. 2*	VB 0714	VB 1014*	VB 1514*
Ordering code:		77100	77300	77500
Research room area	m ²	0.7	1.0	1.5
Dimensions	width	mm	970	1270
	depth	mm	750	750
	height	mm	1400	1400
Climate	Fig. 3	sterile steam humidification		
Temperature range	°C	+5 to +45, without irradiation		
	°C	+10 to +45, with irradiation		
Temperature fluctuation	K	± 0.5 temporally		
Humidity range	%r.h.	40 to 95, without irradiation		
	%r.h.	40 to 85, with irradiation		
Humidity fluctuation	% r.h.	± 3 to 5, temporally		
Dew point range	°C	+4 to +35		
Irradiation module 1	Fig. 5	adjustable in height		
Intensity	μmol/m ² s (klx)	450 (35)		
Distance	mm	200		
Day/night cycle		6 steps		
Lamp type	Fig.7/8	fluorescent tubes		
Air flow module 1	Fig. 6	horizontal with airflow fiber		
Communications module				
Terminal	Fig. 11	C-terminal		
Programming		real time		
Interface		RS 232		
Exterior dimensions	width	mm	1980	2280
	(without door) depth	mm	850	850
	height	mm	1970	1970
Door	width x height	mm	900 x 1400	2x900x1400
Observation window WxH	mm		400 x 700	2x400x700
Nominal voltage		400 V +6/-10%, 3/N, 50 Hz		
Nominal power	kW	3.0	3.3	4.0
Options				
Irradiation module 2	Fig. 4	separate lamp compartment, air-cooled		
Intensity	μmol/m ² s (klx)	700 (50)		
Distance	mm	500		
Day/night cycle		6 steps		
Lamp type	Fig. 9	metal-halide lamps + krypton		
Nominal voltage		400 V, +6/-10 %, 3/N, 50 Hz		
Nominal power	kW	4.5	5.3	7.0
Irradiation module 3	Fig. 4	separate lamp compartment, air-cooled		
Intensity	μmol/m ² s (klx)	1000 (71)		
Distance	mm	500		
Day/night cycle		6 steps		
Lamp type	Fig. 9	metal-halide lamps + krypton		
Nominal voltage		400 V, +6/-10 %, 3/N, 50 Hz		
Nominal power	kW	6.0	7.8	10.5
Height module	Fig. 4	height extension by 500 mm		
Air flow module 2	Fig. 4	vertical, over perforated floor		
Low temperature module	Fig. 4	-10 to +45 °C, without irradiation		
UV irradiation	Fig. 10	on request		
Software TSI	Fig. 12	DOS or OS2		
CO ₂ measuring and dosing unit		on request		

Performance values refer to 25 °C ambient temperature, without test material, without options - We reserve the right of changes in construction resulting from technical progress.

Spectra have been selected according to PAR (Photosynthetic Active Radiation). The combination of different lamps simulate a spectrum similar to that of daylight. The natural simulation of a day is achieved by the discontinuous control of six irradiation intensities. Dawn and dusk are simulated by a separate irradiation intensity which has a high proportion of red light (phytochrome).

The irradiation module 1 - Fig. 5 (450 $\mu\text{mol}/\text{m}^2\text{s}$ at a distance of 200 mm) is equipped with fluorescent lamps which are infinitely variable in height. Adjustable from the outside via an electromotor, the experimental conditions remain constant and undisturbed.

Bio Line offers two other irradiation variations for every size of research room: 700 and 1000 $\mu\text{mol}/\text{m}^2\text{s}$ at a distance of 500 mm (e.g. for C_4 plants). The metal-halide lamps (Fig. 9) are installed in a separate air-cooled lamp compartment. The climate remains constant. Lamps may be easily replaced.

UV radiation (option - Fig. 10) is possible upon request.

Programming of a day/night cycle e.g. short day term or long day term is performed menue-guided in real-time mode with the C-terminal or PC. The microprocessor control monitors the functioning of the system. Faults are displayed on the LCD display of the C-terminal (Fig. 11).

32 systems may be programmed and controlled with the option TSI (Test System Interconnection, Fig. 12). Programmes and research conditions are documented and are graphically displayed. Control limits may be selected.




Fig. 7: Fluorescent lamps - LUMILUX®




Fig. 8: Fluorescent lamps Fluora®




Fig. 9: Metal-halide lamps (option)




Fig. 10: UV-radiation (option)




Fig. 11: C-terminal




Fig. 12: TSI-Software (DOS) (option)

Walk-in Chambers

The concept of the walk-in *Bio Line* chambers is also based on the application orientated modularity. Irradiation intensity, growth height and temperature range as well as the type of air flow can be chosen freely. The standard air flow version is vertical via perforated floor.

Spectra have been selected according to PAR (Photosynthetic Active Radiation). Dawn and dusk are simulated by a separate irradiation intensity which has a high proportion of red light (phytochrome). Natural spectral simulation is achieved by the discontinuous control of six irradiation intensities.

Humidification is performed with sterile vapour. The humidifier requires demineralized water and is maintenance-friendly. A separate, high performance dehumidification system ensures stable climate conditions.

Bio Line does not use materials which contain substances proven to be phytotoxic. The polymeric materials assigned in the research room have been gas-chromatographically analyzed for volatile organic emissions.


The research room is made of insulated elements and is sealed vapour-tight. The insulation material ensures excellent heat insulation which contributes to high constancy of climate conditions and to the reduction of operating costs.

Highly reflecting surfaces create a natural scattering of rays and thus extremely uniform irradiation is achieved over the entire research room.

Bio Line uses the chlorine-free refrigerant. Ozone depletion potential is thus reduced to zero. The insulating material is 75 mm PU-foam which is free from CFC.

Programming of the day/night rhythm is performed menue-guided in real-time mode with the C-terminal (Fig.11) - see control unit of the cabinets. Horizontal air flow (option) with airflow fiber ensures low air velocity and high constancy of climate conditions in the research room.

Bio Line offers a system of options with various functional modules e.g. the irradiation modules and the air flow modules.



Technical Data

Bio Line Walk-in Chambers

Type	VB	4018	8018	12018
Research room, area	m ²	4.0	8.0	12.0
Dimensions	width	mm	2550	2550
	depth	mm	1630	3230
	height	mm	1800	1800

Climate	Fig. 3	sterile steam humification
Temperature range	°C	+5 to +45, without irradiation
	°C	+10 to +45, with irradiation
Temperature fluctuation	K	± 0.5 temporally
Humidity range	%r.h.	40 to 95, without irradiation
	%r.h.	40 bis 85, with irradiation
Humidity fluctuation	% r.h.	± 3 to 5, temporally
Dew point range	°C	+4 to +44

Irradiation module 1	stationary lamp frame	
Intensity	μmol/m ² s (klx)	450 (35)
Distance	mm	200
Day/night cycle		6 steps
Lamp type	Fig.7/8	fluorescent tubes

Air flow module 1	vertical, via perforated floor
-------------------	--------------------------------

Communications module	Fig. 11	C-terminal
Terminal		real time
Programming		RS 232
Interface		

Exterior dimensions	width	mm	2700	2700	2700
	depth	mm	2400	4000	5600
	height	mm	2495	2495	2495

Door	width x height	mm	800 x 1900		
------	----------------	----	------------	--	--

Threshold	mm	395		
-----------	----	-----	--	--

Observation window	WxH	mm	400 x 400		
--------------------	-----	----	-----------	--	--

Separate control unit	WxDxH	mm	800 x 600 x 2100		
-----------------------	-------	----	------------------	--	--

Nominal voltage		400 V, +6/-10 %, 3/N, 50 Hz			
-----------------	--	-----------------------------	--	--	--

Nominal power	kW	8	15	25
---------------	----	---	----	----

Options

Irradiation module 1 / variant	adjustable in height
--------------------------------	----------------------

Irradiation module 2 *	separate lamp compartment, air-cooled
------------------------	---------------------------------------

Intensity	μmol/m ² s (klx)	700 (50)		
-----------	-----------------------------	----------	--	--

Distance	mm	500		
----------	----	-----	--	--

Day/night cycle		6 steps		
-----------------	--	---------	--	--

Lamp type	Fig. 9	metal-halide lamps + krypton		
-----------	--------	------------------------------	--	--

Nominal power	kW	12	20	28
---------------	----	----	----	----

Irradiation module 3 *	separate lamp compartment, air-cooled
------------------------	---------------------------------------

Intensity	μmol/m ² s (klx)	1000 (71)		
-----------	-----------------------------	-----------	--	--

Distance	mm	500		
----------	----	-----	--	--

Day/night cycle		6 steps		
-----------------	--	---------	--	--

Lamp type	Fig. 9	metal-halide lamps + krypton		
-----------	--------	------------------------------	--	--

Nominal power	kW	17	30	41
---------------	----	----	----	----

Height module 1	height extension by 300 mm
-----------------	----------------------------

Height module 2	height extension by 600 mm
-----------------	----------------------------

Air flow module 2	horizontal with airflow fiber
	exterior dimensions in width + 300 mm,
	in height - 300 mm, threshold 95 mm

Low temperature module	-10 to +45 °C, without irradiation
------------------------	------------------------------------

UV irradiation	on request
----------------	------------

CO ₂ measuring and dosing unit	on request
---	------------

shelf systems	on request
---------------	------------

Software TSI	Fig. 12	DOS or OS2
--------------	---------	------------

* The exterior dimension is 3095 mm - Performance values refer to 25 °C ambient temperature, without test material, without options - We reserve the right of changes in constructions resulting from technical progress.

Fig. 13: Scheme - modularity walk-in chambers



Fig. 14: Title - Walk-in chamber at BASF

The cover picture of this brochure shows the agricultural research station at the Ecobiology Division of BASF, Limburgerhof. Data for the registration of pesticides corresponding to EPA and BBA guidelines (with ^{14}C pesticides) is gained here. The results from phytotron and field studies are almost identical. The climate and daylight conditions prevailing in the countries of the five continents with their respective seasonal and diurnal cycles, wind and CO_2 concentration (350 ppm) are simulated. Sunlight is simulated by a xenon reflector-filter system. Programming, control and documenting is performed centrally via computer.

The plants under research include corn, potatoes, rape, rice, soya, tobacco, tomatoes and wheat.

Literature:
"Climate Simulation for Studies on Fate of Pesticides in Plants", 8. IUPAC, Washington D.C., USA 1994



In the exposition chambers at the GSF conditions prevailing in the Bavarian Alps from Spring 1986 to Autumn 1987 were simulated. Forest decline was investigated with cloned spruces. Seasonal and diurnal rhythms were performed in the temperature range of $+26\text{ }^\circ\text{C}$ to $-14\text{ }^\circ\text{C}$ and in the humidity range of 25 % to 98 % relative humidity. The maximum irradiation intensity amounted to $1200\text{ }\mu\text{mol}/\text{m}^2\text{s}$ PAR in a distance of 1000 mm. The chamber was exposed to ozone ($130\text{-}360\text{ mg}/\text{m}^3$). The wind speed reaches to 0.6 m/s. Acid rain was set to a pH value of 3 or 5.6. The area surrounding the roots of the spruces was tempered between $+2\text{ }^\circ\text{C}$ and $+14\text{ }^\circ\text{C}$. Fresh air and discharge air were cleansed with partial filters and activate carbon filters.

The plants under research include beech trees, oak trees, spruces, tobacco, tomatoes and wheat.

Literature:
"Effects of Ozone, Acid Mist and Soil Characteristics on Clonal Norway Spruce - An Introduction to the Joint 14 Month Tree Exposure Experiment in Closed chambers", Environmental Pollution 64 (1990) 189-207.



Fig. 15 and 16: Exposition chambers at the Research Center for Environment and Health, GSF, Munich

In the phytotron of the Max-Planck-Institute for Cultivation Research high growing teosinthes coming from Mexico mature to inflorescence. The spectrum of the metal-halide lamps and an irradiation intensity of $1200\text{ }\mu\text{mol}/\text{m}^2\text{s}$ PAR in a distance of 1000 mm make "nicking" possible. Programming of the day/night cycle is performed in real-time mode. The research room with a height of 3000 mm is equipped with a horizontal air flow system using the airflow fiber. Thanks to its fine-pored structure, the innovative airflow fiber makes extremely low air rates and high constancy of the climate possible. This fiber is micro-biologically and physiologically safe. Natural evaporation behaviour on the stomata is preserved.

Since 1994, this type of airflow fiber has been successfully applied in seven phytotrons at the Max-Planck Institute for Cultivation in Cologne.



Fig. 17: Max-Planck-Institute for Cultivation Research



Fig. 18: Max-Planck-Institute for Cultivation Research



Fig. 19/20: Study Chair of Phytopathology, University of Constance

Experiments are being conducted at the University of Constance with phytopathological fungi e.g. rusts. For inoculation of the plant host, a temperature range of +16 °C to +22 °C and saturated air over a period of 8 hours is necessary. A system of binary nozzles creates 100 % relative humidity in the research room. Decontamination and safety devices are absolutely essential. The research room made of stainless steel may be disinfected with gaseous formaldehyde. The climatic system is cleansed by an additional thermal disinfection system. Filters, Class S (DIN 29 184) with an efficiency of at least 99.97 % clean the discharge air. Underpressure in the research room prevents phytopathogene spores from escaping when the opening of the door.



Fig. 21/22: Department of Genetic Research of Higher Plants, Technical University of Munich

The Study Chair for the Genetic Research on Higher Plants is working on the establishment of barley transformation. The effects of transposable DNA elements and their application in plant genetics as well as the expression of plant proteins in heterologic plant systems are investigated. The rack system with illuminated shelves not only offers adequate irradiation but also a great deal of space. The correct spectrum is achieved by a combination of lamps.

Research cultures are tobacco and Arabidopsis tissue cultures arising from *Agrobacterium tumefaciens*, embryogene callus cultures of barley arising from particle bombardment.



Fig. 23/24: Biophysic Study Group, Technical University of Munich

The Biophysics Study Group of the Technical University of Munich investigated photosynthesis and the water content of agricultural plants e.g. of spring-sown wheat. The challenge was to simulate open land conditions in the laboratory. Realistic field conditions with wind speeds of 3 m/s are simulated in both research rooms. A CO₂ measuring and dosing system is available. Both chambers are connected via TSI software. Programming, control and documentation of "Otto" and "Else" - the code names of the chambers - are performed centrally via PC.

Literature:
 "Simulation of Natural Weather Conditions in Climate Chambers: Possibilities and Limits", Garden Science, 59 (2), Pages 89-94, 1994.

We reserve the right of changes in constructions resulting from technical progress.
Figures partly with options.



Vötsch Industrietechnik GmbH
Umweltsimulation · Wärmetechnik
D-72336 Balingen-Frommern · Beethovenstraße 34
Telefon (07433) 303-0 · Telefax (07433) 303-4112

eMail info@v-it.com
Internet: <http://www.v-it.com>