Orbitl 1900 Heavy Duty Shaker



- Flat or flask clamp platform options
- · Designed for heavy loads and continuous duty
- Safe for use in temperature controlled environments
- · Ideal for liquid cultures and general mixing


Product Information

The Orbit 1900 is a large shaker designed for continuous duty and heavy loads. Its solid construction and counter balanced drive system provide vibration free operation.



Operation of the shaker can be timed or continuous, while the speed is variable. Both time and speed are set digitally. Setting parameters is easy with the analog knob. Values are shown on the large, three digit display.

Three platforms are available to accommodate a variety of different sizes of containers and vials on the Orbit 1900 Shaker. The flat platform is useful for shaking trays, plates, boxes and other containers. A nonslip mat on the platform surface keeps items in place. Flasks are accommodated on the flask platform through the use of clamps (sold separately). Universal platform can be used with various size flasks, bottles and odd shaped vessels. Custom manufactured platforms/clamps are available. Suitable for use in an incubator or cold room, the Orbit 1900 Shaker fits easily inside LabnetIs High Performance 611 Series Incubators.



Speed range	25 - 300rpm
Timer	0 to 99 minutes or continuous
Motion/orbit size	Horizontally circular, 19 mm
Maximum load	22 lb/10 kg
Ambient operating range	+4° to 65°C
Dimensions (W x D x H)	15.3 x 16.8 x 6.2 inches/38.9 x 42.7 x 15.8 cm
Weight	46.2 lb/21 kg
Electrical*	120V~, 60 Hz

* 044	14			
"Utner	voltages	available	upon	request

Ordering Information				
Catalog Number	Item Description	Price		
S2040-1900	Orbit 1900 Digital Heavy Duty Shaker, base only	Request Quote		
S2040-12	Flat platform with nonslip mat (40 x 50 cm)	Request Quote		
S2040-15	Universal platform adjustable securing hardware. Holds a variety of containers.	Request Quote		
S2040-18	Platform for flask clamps (clamps sold separately)	Request Quote		
S2040-85	Clamp for 125 ml Erlenmeyer flask (max. 28)	Request Quote		
S2040-99	Clamp for 250 ml Erlenmeyer flask (max. 20)	Request Quote		
S2040-09	Clamp for 500 ml Erlenmeyer flask (max. 10)	Request Quote		
S2040-01	Clamp for 1000 ml Erlenmeyer flask (max. 6)	Request Quote		
S2040-02	Clamp for 2000 ml Erlenmeyer flask (max. 5)	Request Quote		

Product Specifications